

# Data Quality Issues and Data Review for Low-level Perchlorate Analysis

Nora J. Conlon, Ph.D.

Stephen DiMattei

Region 1, Quality Assurance Unit

# Strategy for evaluating data quality

- ◆ Know what concentration levels are necessary for making decisions
- ◆ Establish Reporting Limits based on calibration standards (not PQLs)
- ◆ Support the Reporting Limits with Method Detection Limit studies
- ◆ Establish meaningful matrix spike samples
- ◆ Review the raw data

# Uses for perchlorate and sources for the environment

## ◆ Solid rocket propellant – missiles

- Matches
- Road flares
- Military Pyrotechnics
- Fireworks

## ◆ Found at military bases

# MMR

## Massachusetts Military Reservation

- ◆ Also known as Camp Edwards and Otis Air Force Base
- ◆ Located over the sole source aquifer on Cape Cod
- ◆ Several plumes affect public and private wells both on and off base
- ◆ Plumes contain mainly TCE, PCE, EDB, or Explosives (RDX) and perchlorate

# MMR (continued)

- ◆ Impact Area is over the highest point of the groundwater table
- ◆ Perchlorate is associated with explosives and propellants
  - Most releases associated with military training and munitions disposal

# Spring 2002

- ◆ Perchlorate sampling and analysis program run by the National Guard Bureau with the ACOE – analyzing to 1 ppb to characterize the site
- ◆ Detections of perchlorate in the wellfield and water supply wells in Bourne, MA at up to 1 ppb
- ◆ Wide spread disbelief that the data were reliable

# QA Steps In – Assessing Data Quality

- ◆ Reviewed a raw data package
  - Chromatograms and calibrations were key
- ◆ Requested MDL studies and SOPs
- ◆ EPA audited the laboratory
- ◆ Proof is in the data
  - The data supported that the lab was able to generate quantifiable data at 1 ppb
  - Samples had little to no interference –  
Low conductivity and low turbidity

# July 30, 2002 QA Memo

- ◆ Jointly issued: RPM and QA
- ◆ Institution of steps to help establish data of sufficient quality to support quantitations at 1 ppb
- ◆ Based on modification of Method 314.0 but does not prescribe methods or instrumentation



# Key Points

- ◆ Establish reporting limit at 1.0 ppb with a low calibration standard
- ◆ Use a 1.0 ppb standard for the MDL study
- ◆ Run a daily MDL check standard at 0.5 ppb (Acceptance criteria 50-150%)

# Key Points - continued

- ◆ Limit method blank contamination to less than the MDL – this requirement applies to pretreatment cartridges as well
- ◆ Use 10 ppb spikes – LCS, MS and MSD (50 ppb spikes are listed in 314.0)

# Current Status

- ◆ Started with a single laboratory capable of meeting these quality control criteria
- ◆ Now have 4 labs that can consistently produce data of known and documented quality at the 1 ppb reporting limit for water
- ◆ Confirmed through split samples and PE samples

# Confidence in the Data

- ◆ The National Guard Bureau (and the Army Corps), Massachusetts and EPA Region 1 have confidence that reliable data is being produced.