

Detection of G-Series and VX Chemical Warfare Agent Degradation Products in Water by Solid Phase Microextraction

¹Robert McKenzie, ²Arthur P. Lee, ^{2*}Philip A. Smith ¹United States Army Center for Health Promotion and Preventive Medicine, Aberdeen MD, ²Uniformed Services University of the Health Sciences, Department of Preventive Medicine and Biometrics, Bethesda MD *Corresponding author

Abstract

- The ability to detect the presence of Chemical Warfare agents in a variety of matrices is important for the protection of both military and civilian populations.
- We looked at the hydrolysis products of the GB, GD, and VX agents using Solid Phase Microextraction (SPME).

Hydrolysis Products

- Isopropyl Methyphosphonate (IMP)
 - GB (Sarin)
 - Isopropyl methylphosphonofluoridate
- Pinacolyl Methylphosphonate (PMP)
 - GD (Soman)
 - Pinacolyl Methylphosphofluoridic Acid
- Ethyl Methylphosphonate (EMP)
 - VX
 - O-ethyl-S-[2-diisopropylaminoethyl] methylphosphonic acid

Instrumentation


■ Gas Chromatograph

- Agilent 5890 Gas Chromatograph
- Merlin Microseal Septum
- Unis 2000 Injection Port

■ Mass Spectrometer

- Agilent 5973 Mass Selective Detector (MSD)
- Inert Source

SPME




Injection Port Liners




SPME Fiber and the Holder (Solid Phase Micro-Extraction)



SPME Fiber Extended



SPME Process



SPME

*SPME Fiber Extended into
the Solution
on Stir Plate
Showing the Fiber
Placement in the Vortex*



GC Conditions

■ GC/MS conditions.

– Chromatographic conditions.

- 15 M x 0.25 mm DB-1701 column with a 0.25 μm Film.
- Injector temperature: 180° C
- 2 minute valve time
- Oven was held at 40° C for two minutes.
- Oven then ramped: 40° C to 180° C at 16° C/minute
- Carrier: Helium with a flow of 2 mLs/minute



Mass Spec Conditions

- Agilent 5973 MSD
- New Agilent Inert Source
 - Improved Sensitivity
 - Improved Peak Shape
- Solvent delay: 2 minutes
- Mass Spectra were collected in SIM mode.

Selected Ions

■ Ethyl Methylphosphonate

- Target Ion 97 Dwell Time 100 μ seconds
- Qualifier Ion 125 Dwell Time 100 μ seconds

■ Isopropyl Methylphosphonate

- Target Ion 97 Dwell Time 100 μ seconds
- Qualifier Ion 79 Dwell Time 100 μ seconds
- Qualifier Ion 125 Dwell Time 100 μ seconds

■ Pinacolyl Methylphosphonate

- Target Ion 123 Dwell Time 100 μ seconds
- Qualifier Ion 124 Dwell Time 100 μ seconds

SPME Fibers

- Tested various fibers:
 - Carbowax (CW)
 - Carboxen/Polydimethylsiloxane (CAR/PDMS)
 - Polydimethylsiloxane/Divinylbenzene (PDMS/DVB)
 - Polyacrylate
 - Polydimethylsiloxane (PDMS)

SPME Test Conditions

- The fibers were tested under four conditions.
 - Water
 - Water with 1 gram of Sodium Chloride
 - Water at pH 4.0
 - Water at pH 4.0 with 1 gram of Sodium Chloride
- The best results were obtained with the addition of Sodium Chloride.

Fibers (with 1 gram of NaCl) vs Mass Spec. Response

Fiber	IMP	EMP	PMP
Carbowax	1056163	770625	1784830
<i>Carboxen/PDMS</i>	5745663	2392218	5335299
<i>PDMS/DVB</i>	16348626	3893553	7222908
<i>Polyacrylate</i>	859319	No Peak	3702256
<i>Polydimethylsiloxane</i>	2644239	859467	5427659

SPME Conditions


- Sample Size: 3 milliliters
- 1 gram sodium chloride added to the sample.
- Extraction time: 30 minutes
- Desorbtion temperature: 180°C
- Desorbtion time: 2 minutes




Retention Times

<i>Compound</i>	<i>Retention Time</i>
<i>Ethyl Methylphosphonate</i>	5.19 minutes
<i>Isopropyl Methylphosphonate</i>	5.58 minutes
<i>Pinacolyl Methylphosphonate</i>	8.39 minutes


EMP Chromatogram




IMP Chromatogram




PMP Chromatogram




EMP Standard Curve



IMP Standard Plot




PMP Standard Plot




EMP Timed Exposure Study

**Ethyl Methylphosphonate Fiber
Exposure Times**




IMP Timed Exposure Study

**Isopropyl Methylphosphonate Fiber
Exposure Times**



PMP Timed Exposure Study

**Pinacolyl Methylphosphonate Fiber
Exposure Times**



Method Detection Limits

- Method Detection Limits (MDLs).
 - Used the EPA Method to determine the MDLs.
 - 40 CFR 136, Appendix B
 - 9 replicates were used,

<i>Compound</i>	<i>MDL</i>	<i>Standard Deviation</i>	<i>Average Recovery</i>
<i>Ethyl Methylphosphonate</i>	29 µg/L	10 µg/L	104 %
<i>Isopropyl Methylphosphonate</i>	0.14 µg/L	0.048 µg/L	99 %
<i>Pinacolyl Methylphosphonate</i>	4.50 µg/L	1.5 µg/L	84 %

Conclusion

- The calibration curves are linear with an r^2 value of at least 0.995.
- The analytes can be seen in the low $\mu\text{g/L}$ range.
- The Analysis is very technique dependent.
 - Reproducibility of the placement of the fiber in the sample vortex is critical.
 - Timing of the fiber exposure has to be reproducible.